

ProMTag reversible protein labeling agent

Catalog No:	T-6000
Molecular weight:	529.55 g/mol
Concentration:	30 mg/mL, 56.65 mM
Shipped in:	90% acetonitrile, 0.1% trifluoroacetic acid
Appearance:	Pink liquid
Storage:	Upon receipt store at -20 °C
Shelf life:	Three years at -20 °C

About:

ProMTag is a pH reversible protein labeling agent. At pH 8.0 or higher, ProMTag reacts with amines on the surface of proteins to covalently link an affinity tag to the surface. At pH 3.0 or less, this covalent link is broken, returning the protein amines to their original, unmodified state. This reversible tag can be used for fast affinity purification without the need for permanent covalent or genetic tags that could interfere with the native mass or isoelectric points. Pair this tag with the ProMTag Capture Resin (Catalog No. CSB-8000) for optimal results. To also digest proteins into peptides, we also recommend using mP-trypsin (mPT-9000), which is stabilized and ProMTag labeled to ensure complete digestion and removal of trypsin after digestion.

Questions about compatibility, optimization, or troubleshooting? Email <u>info@impactproteomics.com</u> for help.



General ProMTag protocol:

- Make sure your protein sample starts in a buffer pH 8 or greater that does not contain primary amines. Our favorite buffers are 100 mM HEPES, pH 8.0 and 100 mM Sodium Phosphate, pH 8.0.
- Dilute your starting sample to 1 mg/mL. For every 10 μg of protein, add 25 μg of ProMTag. Vortex briefly to mix
 - a. Ex: If you are labeling 100 μg of protein, you need to add 250 μg (8.3 $\mu L)$ of ProMTag
- 3. Incubate for 30 minutes to label your protein sample.
- 4. Add sample to the appropriate amount of ProMTag Capture Resin to capture labeled proteins. Incubate for 30 minutes with end over end rotation to keep the beads suspended.
 - a. For every 10 μg of labeled protein, you need 20 μL of ProMTag Capture Resin.
 - b. If you are using a buffer that has SDS, you need at least 10 μL of ProMTag Capture Resin for every 10 μL of 1% SDS. The total amount of SDS is what is important here so 10 μL of resin will also tolerate 20 μL of 0.5% SDS or 5 μL of 2% SDS.
- 5. If you are washing away contaminants or doing a buffer exchange, wash the resin at least three times with a wash buffer that is at pH 8.0 or greater. Use a volume of wash buffer that is at least five times the volume of capture resin if possible. Ex: For 40 μL of beads, wash with 200 μL of wash buffer.
 - a. If washing away nucleic acids, we recommend adding a wash step that uses 250 mM sodium chloride.
 - b. If washing away SDS, we recommend adding a wash step that uses a wash buffer that contains 50-75% acetonitrile
 - c. For the final wash, wash using just purified, deionized water
- 6. To elute proteins from the ProMTag Capture Resin, add a buffer with a pH of 3.0 or less at a volume equal to the volume of ProMTag Capture Resin



used. Our favorite buffer uses 100 mM formic acid. Incubate for 15 minutes to release the protein from the resin

- a. Ex: If you used 40 μL of ProMTag Capture Resin, you need 40 μL of elution buffer.
- b. If you would like to digest proteins into peptides, do not centrifuge your sample to elute, only add the elution buffer and incubate for 15 minutes. Then continue on to the next step.
- c. If you are not digesting proteins, repeat the elution step once to ensure all protein is eluted.
- If you would like to digest proteins into peptides using mP-trypsin, add 10 μg of mP-trypsin to the sample for every 10 μg of protein that needs to be digested. Incubate for 1 h to digest.
 - a. During this time the mP-trypsin will also covalently attach irreversibly to the ProMTag Capture Resin to ensure a clean final sample
- 8. After digestion, spin to elute the final sample. Add the same volume of elution buffer and elute one more time to ensure capture of the entire protein sample.