

Reversible chemistry for universal protein extraction and cleanup of whole proteins and peptides for fast, high yield proteomics sample preparation

Stephanie Biedka¹; Jonathan Minden^{1, 2}; Nolan Frey²; Amber Lucas¹


¹Impact Proteomics, Pittsburgh, PA; ²Carnegie Mellon University, Pittsburgh, PA

Abstract


High quality protein extraction and cleanup during sample preparation is critical for achieving the coverage, yields, and reproducibility required for a successful proteomics experiment. Current technologies for proteomics sample preparation generally rely on precipitation or filtration-based technologies, both of which can suffer from sample loss, long processing times, incompatibility with certain buffers that aid in extraction, and limitations for automation. Recently, new sample preparation technologies have made improvements to filter- and precipitation-based sample preparation and have come closer to creating efficient protein sample preparation workflows. However, these methods still suffer from many of the aforementioned issues, which can only truly be addressed by taking a novel approach to protein sample preparation

Biochemical tags have been used for decades for preparation of single protein samples, with affinity chromatography being a fast and simple way to separate a target protein from other cellular material. Biochemical tags have not previously been used for proteomic sample preparation in a similar manner due to numerous technological challenges, such as the tags altering protein size and charge, making mass spectrometry analysis complicated or impossible. We have created a covalent, reversible chemical tag, which we refer to as ProMTag, that is used for tagging, capture, and cleanup of whole proteome samples using a complimentary resin that allows for extensive washing without sample loss. This reversible tagging method allows the convenience of affinity purification without permanent modification of the protein sample, ensuring fast, high yield, and easy to perform sample preparation that yields highly reproducible results and requires less than 30 minutes of hands on time. With this novel approach to sample preparation, this technology can be used to cleanup samples for single protein projects, whole protein gel-based proteomics workflows, or peptide mass spectrometry-based proteomics workflows without compromising yields or coverage.


Universal Protein Extraction and Cleanup Kit (UPECK) workflow




ProMTag reversibly reacts with proteins




ProMTag labeled protein can be captured on ProMTag capture beads




Proteins eluted from ProMTag capture beads in 15 minutes




Whole protein proteome cleanup using UPECK



Protein digestion and Trypsin removal using mP-12 Trypsin



Yields and mass spectrometry of UPECK prepared samples



Yields for peptides generated using UPECK were assessed using the Pierce™ Quantitative Fluorometric Peptide Assay in various popular buffer systems. All buffer systems averaged yields around 85% for 20 µg of starting sample.

Yeast proteomes were cleaned up using four different techniques and were sent for mass spectrometry analysis. UPECK showed more identifications of proteins overall and more unique protein identifications. Additionally UPECK showed a lower false discovery rate and higher protein coverage than the other techniques. Note: Because of Covid-19, we were only able to perform this once. This needs to be repeated for statistical significance.

Acknowledgements

We would like to acknowledge and thank the Vincent Coates Foundation Mass Spectrometry Laboratory at Stanford University Mass Spectrometry for perform mass spectrometry analysis for us. This work was supported in part by NIH P30 CA124435 utilizing the Stanford Cancer Institute Proteomics/Mass Spectrometry Shared Resource. This work was also supported by NSF SBIR Phase I award 1843332.

