

# Fast, High-Yield, and Universal Proteomics Sample Preparation Using the Reversible Protein Tag, ProMTag

Svitlana Yablonska<sup>1</sup>, Stephanie Biedka<sup>1</sup>, Charity L. Washam<sup>2</sup>, Xi Peng<sup>3</sup>, Stephanie Byrum<sup>2</sup>, Kevin Xiao<sup>3</sup>, Jonathan S. Minden<sup>1</sup>

<sup>1</sup>Impact Proteomics, Pittsburgh, PA, <sup>2</sup>UAMS, Little Rock, AR, <sup>3</sup>Center for Clinical Mass Spectrometry, AHN Cancer Institute, Pittsburgh, PA



## Abstract

High-quality protein extraction and cleanup during sample preparation are critical for achieving the coverage, yields, and reproducibility required for a successful proteomics experiment. Many current technologies for proteomics sample preparation rely on precipitation or filtration-based technologies, which suffer from sample loss, long processing times, incompatibility with certain lysis buffers, and limitations for automation.


At Impact Proteomics, our Universal Protein Extraction and Cleanup Kits (UPECK) have overcome these limitations. Our novel protein sample preparation method utilizes a protein-specific tag (ProMTag) where one end of the tag forms a reversible, covalent link to proteins, and the other end irreversibly anchors the protein to a capture resin using a fast-acting, click chemistry reaction. The ProMTag reaction is not affected by high detergent, denaturant, or salt conditions accommodating harsh cell lysis conditions. Once ProMTagged proteins are covalently linked to the capture resin, contaminating salts, detergents, and other molecules that interfere with MS analysis are removed by rapid washing steps. Following washing, the ProMTag is reversed, releasing the protein in its original, unmodified form. To yield peptides ready for MS analysis, modified MT-Trypsin, which also binds to the capture resin, is added and digestion is complete in one hour.

The UPECK workflow ensures fast, high-yield, and streamlined sample preparation that yields highly reproducible results and requires about four hours to complete with less than an hour of hands-on time.


Here, we demonstrate the compatibility of the UPECK Peptide workflow with various sample types, including tissues from brain to muscle, and bodily fluids, all using a universal lysis and solubilization buffer.

## ProMTag UPECK-Pep ensures streamlined sample preparation


- A reversible, covalent link is formed by one end of the ProMTag, and an irreversible anchor is formed by the other end via a fast-acting click chemistry reaction.
- Rapid washing steps are used to remove contaminating salts, detergents, and other molecules that would interfere with MS analysis.
- The ProMTag is reversed after washing, allowing the protein to return to its original state.
- A modified MT-Trypsin, which also binds to the capture resin, is added to yield peptides ready for MS analysis.




## Peptides were prepared in triplicate via ProMTag UPECK-Pep from a variety of samples



## Triplicate samples prepared with ProMTag UPECK-Pep showed high reproducibility




## Analysis confirms high consistency of sample replicates




## The ProMTag UPECK-Pep workflow is compatible with a wide variety of sample types

### PCA of MS data shows high reproducibility and similarity between samples



### Pearson's correlation analysis confirms a high correlation between sample replicates and tissue origins



## Acknowledgements

Funding for this work was provided by the NSF SBIR Phase II Award 2036199 (S.Y., S.B., J.S.M.), NIH grant P20GM121293 and the UAMS Winthrop P. Rockefeller Cancer Institute (C.L.W., S.B.).