

Abstract


Autoimmune diseases affect >20 million people in the US today. Currently, disease-specific autoantibodies are thought to be the best biomarkers for diagnosis. Conventional immunoprecipitation methods have been used to identify autoantigens from the most common autoimmune diseases. However, these diseases account for only 6.5 million of the 20 million patients suffering from autoimmune diseases, leaving many without diagnoses until irreversible damage occurs. The remaining 13.5 million patients have >80 autoimmune disorders without well characterized autoantibodies. The state-of-the-art diagnostic test of these remaining diseases relies on gel electrophoresis of immunoprecipitated radiolabeled proteins, which cannot be identified by MS due to safety issues and the overwhelming presence of immunoglobulins.

We have created an immunoprecipitation method that uses serum from patients with any autoimmune disorder to identify patient-specific autoantigen proteins. This method uses a reversible click chemistry tag, called ProMTag. One end of the ProMTag forms a reversible, covalent bond with protein by coupling to lysines and amino termini. The other end of the ProMTag can form an irreversible, covalent bond with a solid bead support via a click chemistry, methyltetrazine-TCO, pairing. In this study, the proteins of cell lysates that contain potential autoantigens were labeled with ProMTag. The ProMTagged-proteins were exposed to patient antibodies bound to Protein A beads, thus capturing the ProMTagged autoantigens. All proteins were released from the Protein A beads, including ProMTagged- autoantigens and untagged-antibodies. The ProMTagged-autoantigens were subsequently coupled to TCO beads, and the untagged-antibodies were washed away. The linkage between the ProMTag and autoantigens was then reversed, yielding autoantigen proteins with greatly reduced antibody contamination ready for MS analysis. MS analysis successfully identified autoantigens from patient serums with scleroderma or myositis/anti-synthetase syndrome. This autoimmune biomarker discovery method can accelerate sample testing for known autoantigens and facilitate rapid discovery of novel autoantigens for both diagnostic and predictive biomarkers.

ProMTag workflow key components

ProMTag Protein tagging


Reversible Protein Binding Reagent
Carboxyethyl-dialkyl-maleic Anhydride (CDM)




Irreversible Bead Binding Reagent
Methyltetrazine (MT)

ProMTag Capture TCO-bead


Trans-cyclooctene (TCO)



ProMTag immunoprecipitation method for purification of autoantigens



Gel-based validation of the ProMTag immunoprecipitation workflow




Control immunoprecipitation of K562 lysate with normal serum

Normal Patient Serum Control					
Description	Log Prob.	Total Intensity	# of spectra	# of unique peptides	Coverage %
Albumin	447.63	1.72E+09	197	69	78.82
Complement C3	443.94	9.78E+08	211	85	54.66
Vimentin	338.95	1.05E+09	135	50	74.03
Apolipoprotein B-100	224.01	1.62E+08	53	46	14.42
CD5 antigen-like	164.28	4.40E+08	72	22	72.05
Alpha-2-macroglobulin	84.98	8.10E+07	24	17	25.44
Polyadenylate-binding protein 1	81.26	1.39E+08	27	15	30.66
Tubulin beta chain	69.21	4.70E+07	13	8	32.43
Actin, cytoplasmic 2	56.92	1.37E+08	22	14	52.00
Complement C4-B	42.44	1.41E+07	5	4	3.10

Mass spectrometry results of negative control identify background proteins that are in the sample after the ProMTag immunoprecipitation workflow.

Immunoprecipitation of K562 lysate with four characterized patient sera



Autoantibodies from patient sera with known autoantigens were immobilized on Protein A beads and used for immunoprecipitation of a ProMTag labeled protein lysate from K562 cells. P1 and P2 had myositis/anti-synthetase syndrome; P2 and P4 had scleroderma (anti-Topoisomerase I). Immunoprecipitated antigens were then sent for mass spectrometry identification. In all four tests, the correct autoantigen was identified in the mass spectrometry result. Notice that P4 had an additional, novel autoantigen.

Conclusions

The ProMTag immunoprecipitation method allowed for mass spectrometry analysis of autoantigens immediately after immunoprecipitation, something that previously could not be done because of antibody contamination. Additionally, we correctly identified all known antigens in clinically relevant patient samples.

Acknowledgements

This work was supported by NSF SBIR Phase II Award 2036199. This work was also supported in part by NIH grant P20GM121293 and the UAMS Winthrop P. Rockefeller Cancer Institute and by NIH P30 CA124435 utilizing the Stanford Cancer Institute Proteomics/Mass Spectrometry Shared Resource at the Vincent Coates Foundation Mass Spectrometry Laboratory, Stanford University Mass Spectrometry (RRID:SCR_017801). We thank Ryan Leib and the Stanford University Mass Spectrometry Facility for their help with mass spectrometry analysis, data analysis, and statistical analysis.